

NEW: LABELING of Functional Tests and Analyses

As your reliable and innovative partner for diagnostic services, we are committed to constantly expanding our offerings for you and improving our processes for you and your patients. Sometimes, it's the little things that can make a difference in the daily laboratory routine. Today, it is a sticker.

With this, you can now label your samples for all functional tests – preferably prominently visible on the shipping bag.

Your Benefits

- ✓ Your highly sensitive samples reach the right laboratories on a daily basis
- ✓ Direct processing without delays
- ✓ Improved sample stability through faster processing



IMPORTANT! Please use the sticker exclusively for the labelling of the following tests:

- **LTTs** (C650, D1010, D1020, D1025, D1060, D1070, D1080, D1090, K115)
- **ELISpots** (K120, K121, K221, K265, K236, K335, K346A, K242, K260, K271, K505, K620)
- **NK-Cell-Test** (D500, D510)
- **BHI**
- **+ Extended Parameters** (E328, E335, E336, E337, E338, E339)
- **Mitochondrial Activity** (E330)
- **Cytokine Status** (D845, D850)
- **TNF- α -Inhibition Test** (D865, D8650)

The existing labelling of the shipping bags for the EDIM Apo10/TKTL1 with a red sticker will remain unchanged. The sticker can now be ordered via our material order form: "Sticker for functional tests", item designation: **IMMUAUF**.

We thank you very much in advance and send you
warm Regards from Limburg

Your biovis Team

Biovis Diagnostik MVZ GmbH
Brüsseler Str. 18
65552 Limburg-Eschhofen
Tel. +49 (0)6431 21248 0
Fax. +49 (0)6431 21248 66

www.biovis.de
Bank: Kreissparkasse Limburg
IBAN: DE12345678945
BIC: ABCDEFG
USt.-IdNr.: DE123456789